06.75.19

CNEA AC 6/75

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ANEXO 1

Contamination Monitor H 1359 A

optionally with gamma compensation


The H 1359 contamination monitor is designed to give you still greater safety in your isotope laboratory. With it, contamination on the hands, footwear, clothing, work surfaces, floors, etc., is located quickly and accurately. The lightweight, large-area HGZ I and III counters and the general handiness of the unit itself make this monitor equally suitable


for both stationary and portable use. The various accessories and detectors available for the instrument allow it to be used in radiation protection as floor monitor, clothing monitor, hand-and-foot monitor and personnel monitor, and also as a neutron ratemeter.

Principal Features

- Simple and functional in use
- Lightweight and compact
- Suitable for a wide range of duties
- Built-in counting-gas supply
- Built-in lead storage batteries for operation independently of the mains
- Mains/buffer-battery operation with electronic-controlled charging unit
- Electronics with integrated circuits, exchangeable plug-in cards
- Six linear measuring ranges and three time constants
- Acoustic signalling of individual pulses
- Warning threshold with optical and acoustic indication
- Built-in Ra control preparation for alpha and beta test
- Recorder connection for continuous recording
- Optionally with differential circuit for gamma compensation (H 1359 K)

HGZ I

HGZ III

Attachable Detectors and their Applications

- Large-area Counter HGZ I, counting surface 180 cm² handy detector for measurement of smallest alpha and beta activity
- Tritium Counter HGZ III, counting surface 11 cm² very light and handy large-area counter, specially suitable for measurement of tritium contamination
- Large-area Counter HGZ II, counting surface 600 cm² for measurement of smallest alpha and beta activity; suitable for stationary use, e.g., as hand-and-foot monitor, floor monitor, clothing monitor, air monitor, water monitor, etc.
- Large-area Counter HGZ I/K and HGZ II/K with Screen Counter, counting surface 180 and 600 cm² when used in conjunction with the H 1359 K for gamma compensation by subtraction, especially in monitored areas with increased environmental radiation
- Neutron rem Counter 5210 A for dose-rate measurement of neutrons from the thermal range to 17 MeV with an approximately ideal rem pattern in accordance with ICRP requirements
- G.M. gamma counter tubes for measurement of gamma and X-ray (roentgen) radiation
- G.M. end-window counter tubes for location of point-pattern contamination
- BF₃ counter tubes for measurement of thermal neutrons
- Proton-recoil counters for measurement of fast neutrons

Description

In the H 1359 Contamination Monitor you have a handy and simple-to-use radiation meter for detection of the slightest specific alpha and beta activity. Radioactive contamination on hands, footwear, clothing, work surfaces, floors, etc., are located rapidly and with certainty. In this, the lightweight, large-area counters and the compactness of the instrument itself (which can be operated either from the mains or from the built-in lead storage batteries), are an advantage.

For measuring contamination you have a choice between the large-area counter HGZ I or the compensated version HGZ I/K, with an effective counting surface of 180 cm² or — specially for tritium contamination — the HGZ III counter with an effective counting surface of 11 cm². If the instrument is used as a stationary monitor you can also use the HGZ II or HGZ II/K counters, with an effective counting surface of 600 cm². These cast aluminium counters are strong and rugged, light in weight and easily decontaminated.

In addition to these large-area counters, you can also use the monitor with the rem counter* for neutron measurement (thermal to 17 MeV) and other proportional counters such as BF₃ and proton-recoil counters, as well as our highly successful G.M. counter-tube probes, including end-window counter-tube probes for exact location of point-

pattern contamination, and a series of gamma probes.


The counting-gas for the HGZ I/K and III is supplied from a gas cylinder in the instrument. The cylinder charge is sufficient for about 25 measurements. The counter is flushed by pressing a button on the front panel of the monitor; this opens the gas cylinder valve, and the butane-propane gas mixture flows through a hose into the counter. The HGZ I counter is provided with an automatically-opening outlet valve and can be operated for several hours on a single gas filling.

This uncomplicated method of gas supply allowed us to design a light and handy monitor and detector. If the monitor does not have to be moved from place to place it can be supplied with the gas from larger steel cylinders.

A practical advantage is that you can check the readiness of the monitor and detector at any time for both alpha and beta radiation by means of the built-in Ra control preparation.


When not in use the standard HGZ I counter can be inserted in its protective holder on the instrument.

* L. Ö. Andersson and J. Braun, A Neutron Rem Counter, „Nukleonik“, Vol. 6, No. 5 (Sept. 64) pp 237–241

HGZ II/K

Neutron rem Counter 5210 A

The large counting surface and the very thin window sheeting (not more than 0.7 mg/cm^2), ensure rapid detection of alpha and beta contamination, the values being shown on a large, easy-to-read scale. Detection is made even easier by the built-in loudspeaker which can be distinctly heard even in a noisy environment.

The monitor has a steplessly-adjustable warning threshold, with a red indicator lamp and an acoustic signal which are triggered as soon as the set alarm level is exceeded.

The high voltage for the large-area counters can be steplessly adjusted by potentiometer between 1,000 and 4,000 V, and separately for the alpha and beta plateaus. A selector switch is provided for the various modes of operation — alpha or alpha-beta measurement or G.M. counter tube. The high voltage is electronically stabilized and can be read on the instrument scale. A separate input with a stabilized voltage of 560 V is available for operation of our counter-tube probes.

In the case of operation from the mains, supply is from an electronically controlled constant-voltage power pack with automatic current limitation. The monitor can be switched instantly to mains operation if the lead storage batteries go flat, and at the same time the batteries will be recharged. The constant-voltage power pack protects the batteries against overcharging. The monitor is therefore excellently suited for continuous buffer-battery operation.

Use as hand-and-foot, clothing, and floor monitor

With the various detectors and accessories available, you can use both the H 1359 monitor and the compensated H 1359 K version for a wide variety of radiation measurement work.

A fully effective, yet simple and economical hand-and-foot monitor station requires only the addition of a bracket for mounting the instrument and the HGZ II hand counter, and a platform with grating and the HGZ II foot counter.

For use as a floor monitor, a mounting frame for two counters — HGZ II or HGZ II/K — is available. The distance between the counters and the floor is adjustable, which is an important feature when measuring alpha contamination.

For the monitoring of clothing, large-area counters arranged one above the other give unfailing detection of contamination.

H 1359 K with differential circuit for gamma compensation

In cases where measurement of minimal activity is frustrated by a high background — which occurs particularly in places with above-normal environmental radiation — the use of our compensated counters HGZ I/K or II/K is recommended. These large-area proportional counter tubes consist of two counters one above the other. The upper one has a thin window and is therefore sensitive to alpha, beta and gamma radiation, while the lower screen counter only measures gamma radiation.


With these large-area counters and the H 1359 K monitor — which has two differential ratemeter circuits — that part of the background detected by the screen counter can be subtracted from the pulse rate of the measuring counter. This effectively compensates for the interference radiation and substantially improves the beta detection performance. The compensation can be adjusted by potentiometer.

Electronic system


The electronic system makes extensive use of integrated circuits, which ensures outstandingly good stability. All electronic components are arranged on exchangeable and easily accessible plug-in cards. The pulses from the proportional counter pass to the integrated circuit of the input discriminator, while another integrated circuit serves as monovibrator for the generation of unit pulses. The integrator consists of an operational amplifier, likewise the comparator for the alarm threshold.

A gated transistor oscillator stage provides the high voltage, a certain fraction of which passes from a voltage divider to one input of a monolithic comparator. Connected to the second input is a reference voltage (alpha-beta plateau), which is derived from a reference diode and is separately adjustable from two potentiometers. As soon as the high voltage deviates by about 1 per cent from the set point, the oscillator stage will be immediately gated to counteract the deviation.

The lead storage batteries are charged from an adjustable constant-voltage power pack with adjustable current limitation and consisting of variable-gain amplifier, limiter, reference element and power transistor in an integrated circuit.

H 1359 with HGZ I

Simple hand-and-foot monitor

Specification

Ranges: 0–20; 0–2x10²; 0–2x10³; 0–2x10⁴; 0–2x10⁵; 0–2x10⁶ pulses per min.

Time constants: 1, 10 and 100 seconds

Input sensitivity: approx. 2mV, negative pulse

Warning: alarm threshold steplessly adjustable, indicated by red lamp and audible signal

High voltage: alpha plateau continuously adjustable from 1 to 4 kV; beta plateau continuously adjustable from 1 to 4 kV; detent-type switch to change over from alpha to alpha-beta operation

G.M. range approx. 560 V

High-voltage indication: 1 to 4 kV on meter scale

Power supply: 1). from mains 220 V, 50 Hz; buffer-battery operation; batteries recharged simultaneously; constant-voltage power pack with automatic current limitation prevents overcharging. Approx. 12 hours charging time for 80 per cent of rated capacity

2). from built-in lead storage batteries, about 12 hours on continuous operation

Battery check: on measuring instrument; green range indicates battery in working order

Function check: by built-in alpha-beta-gamma source, approx. 50 nC Ra.

Recorder output: 100 μ A, 0–1.5 kOhm, for "Miniscript" strip recorder or "Elmes" circular-chart recorder

Detector connections: two sockets (DIN 44424) for proportional counters; Cannon socket XLR-5-31 for G.M. counter-tube probes

Dimensions: height 110 mm, width 300 mm; depth 220 mm

Weight: 4.4 kg

Attachable Detectors

Large-area Counters HGZ I and I/K

Effective window surface: 180 cm²

Window weight: not more than 0.7 mg/cm²

Operating voltage: alpha approx. 2.8 kV; beta approx. 3.4 kV (butane-propane)

Length of plateau: alpha approx. 250 V; beta more than 300 V; plateau slope 2 %/100 V

Background: alpha approx. 1 pulse/min.; beta approx. 500 pulses/min.

Gas supply: built-in butane-propane gas cylinder (160 cm³), pushbutton-operated; automatically-opening outlet valve on the counter

Operating pressure: atmospheric (760 Torr.)

Connecting cable: incl. gas connection, length 1.1 metre

Dimensions: base surface 152 x 227 mm; height 23 mm, or 65 mm with handle

Weight: 760 g

Tritium Counter HGZ III

Window thickness: not more than 0.5 μ m

Unit weight: max 70 μ g/cm²

Effective window surface: 11 cm²

Dimensions: base surface 40 x 150 mm; height with handle 50 mm

Operating voltage: 3.6 to 3.8 kV

Weight: 300 g incl. cable and plug

Window sheeting: mounted on an easily exchangeable, reusable polyester panel

Gas supply: from butane-propane gas cylinder (160 cm³) in the monitor, pushbutton-operated

Large-area Counters HGZ II and II/K

Effective window surface: 600 cm²

Window weight: not more than 0.7 mg/cm²

Operating voltage: alpha approx. 1.6 kV; beta approx. 2.1 kV (argon-methane)

Length of plateau: alpha approx. 250 V; beta more than 300 V; plateau slope 3 %/100 V

Background: alpha less than 1 pulse/min.; beta approx. 1,400 pulses/min.

Dimensions: base surface 225 x 390 mm; height 38 mm

Weight: 3.5 kg

Neutron rem Counter

rem-proportional from 0.025 eV to 17 MeV, sensitivity (6 pulses/sec.) / (2.5 mrem/hour)

Proton-recoil Counter

FN2/6a, 20th Century, neutron sensitivity (1.8 pulse/sec.) / MPD

G.M. Counter Tube Probes

Gamma probe IIa	VALVO 18503 with filter
Alpha-beta probe III	VALVO 18506
Gamma probe VI	VALVO 18520
Gamma probe VII	VALVO 18509

Also available are:

hand-foot-clothing monitors, monitors for linen, etc step-filter air-monitor equipments, plutonium aerosol monitors, portable neutron ratemeters with rem counter, multi-channel monitoring systems for rack mounting, semiconductor detectors, He³ counter tubes, AEC/NIM electronics


We also manufacture high-frequency, ultrasonic and heat-impulse welding apparatus for sealing radioactive samples, waste, etc. Please send for descriptive literature.

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Precision-high voltage-unit, type HSP 2

Application

The precision-high voltage-device is intended for the operation of radiation detectors of our proportional counters, and also for Geiger-Müller counters. Its high accuracy and constancy guarantee sure and accurate measurements by simple operating.

Construction

The device is built in a $\frac{1}{4}$ -unit in 19"-system. The unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. All parts are easily accessible and interchangeable.

The controls on the panel are: 1 ten-turn potentiometer and 1 trimming potentiometer for switchable high voltages. Also on the panel there is the indication meter for high voltage.


Performance

An AC generator feeds its voltage in a regulating unit. The regulating unit is controlled by the forward-backward control. The AC voltage supplied from the regulating unit is fed in an amplifier. The output voltage of the amplifier is transformed up. This high voltage passes a voltage multiplier that rectifies it simultaneously. A part of the rectified high voltage is compared with a reference voltage. The DC voltage supplied from this stage controls the regulating unit.


Technical data

(subject to change without notice)

High voltage, adjustable continuously	0...5000 V, positive against ground
Ripple voltage	below 0,05 V at 5000 V
High voltage, linearly adjustable by ten-turn potentiometer	
Deviation from linearity at 500 V	2 ‰
Deviation from linearity above 500 V	1 ‰
2 high voltages, each adjustable with	α , $(\alpha + \beta)$ - switch
1. high voltage adjustable with	ten-turn potentiometer
2. high voltage adjustable with	ten-turn trimming potentiometer
Long time stability	1 ‰
Backward control by voltage reference	
Control meter indication, accuracy	2,5 ‰
Built-in, regulated power supply, intern operating voltage	± 15 V
Line voltage	220 V, 50 cps
Power input	10 VA
Wire fuse	0,2 A
Size	$\frac{1}{4}$ -unit, 19"-system
Height	178 mm
Weight	about 3000 grammes

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
Difference ratemeter unit, type DRM 2

Application

The difference ratemeter unit is intended for forming the difference of two pulse rates. This difference is necessary using continuously monitoring of water and also monitoring of air in rooms with specific measurement of tritium and, higher-energy gases. Also the unit may be used, if extensive independence of varying background and fully-automatic operating or remote control are required.

Construction

The device is built in a $\frac{1}{2}$ -unit in 19"-system. The unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. All parts are easily accessible and interchangeable. The controls on the panel are: electric zero shift, measuring selector switch, compensation fine adjusting knob, switch for connection to automatic range switching, 8 push buttons for measuring ranges, 6 push buttons for time constants, two adjustable warning thresholds on the measuring contactor.

Performance

(see block diagram)

Both pulses supplied from the measuring counter and pulses supplied from the compensation counter are passed on to the measuring selector switch. Before the pulses pass a binary chain for coarse adjustment of the compensation. This binary chain can reduce the compensation pulse rate in large limits. The measuring selector switch permits, beside of the adjustment a full compensation, to make visible the measuring rate and the compensation rate. Both the ratemeter measuring channel and the ratemeter compensation channel possess ranges switched by relays. Hence they are switchable by hand, with remote control or with automatic range switching. A compensation fine adjustable by a tenturn potentiometer completes so the coarse adjustment, that the compensation, i. e. the adjustment of the zero voltage, may be made on the capacitor C, that determines the time constants. The capacitor voltage is amplified by an indication-record amplifier and fed in the record output and measuring contactor. The latter indicates the measuring value and has minimum and maximum contacts for the adjustment of thresholds.

Technical data

(subject to change without notice)

Input pulses from	$\frac{1}{4}$ -unit V 22, AV 2, GV 2 or AKV 2
8 linear measuring ranges	100...300, 1000...3000, 10000...30000, 100000...300000 pulses per minute
Accuracy of indication	2 %
Deviation from linearity within a range	1 %
6 time constants, adjustable by push buttons	2, 4, 16, 32, 64, 128 sec.
Indicator with 2 independent contacts for adjustment of a maximum and minimum warning threshold	measuring contactor
Measuring rate M and compensation rate K	each adjustable
Compensation factor	
Coarse adjustment $\begin{matrix} K \\ M \end{matrix}$	1 : 1 to 16 : 1
Fine adjustment $\begin{matrix} K \\ M \end{matrix}$	by ten-turn potentiometer
Modes of operation: Measuring ranges to select by hand, remote control, and automatic range switching	
Electric zero shift, dis- connectable, fixed	+ 20 %
Record output	max. 10 V, max. 0,01 A against ground
Power supply from	$\frac{1}{4}$ -unit STW 2
Size	$\frac{1}{2}$ -unit, 19"-system
Height	178 mm
Weight	about 4300 grammes

Block diagramm DRM 2

(from top left to bottom right)

measuring selector switch — ratemeter measuring channel
— time constant C — indication-record amplifier — range
switching — measuring contactor — binary chain — rate-
meter compensation channel — compensation, fine — mini-
mum-maximum threshold —

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
Power supply-warning-unit, type STW 2

The device is designed in a $\frac{1}{4}$ -unit of 19"-system and intended for the power supply of all units in our production line. This unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. The device supplies 3 short-circuit-proof supply voltages. It contains a tone oscillator that provides a voltage of higher or lower frequency, if connected from outside. By a built-in loudspeaker the pulses from 2 channels (e. g. measuring channel or compensation channel) or a higher and lower warning tone are made audible by means of a selector switch. The volume is adjustable.


Technical data

(subject to change without notice)

Line voltage	220 V, 50 cps
Power input	about 40 VA
Wire fuse (on the rear of unit)	0,3 A (time-lag)
Supply of 3 regulated and short-circuit-proof voltages	+ 15 V (0,5 A) + 5 V (0,2 A) - 15 V (0,5 A) against zero potential
Tone oscillator, to be connected from outside	high tone about 2000 cps low tone about 400 cps
Loudspeaker with adjustable volume	
Tumbler switch with 3 positions for listening of pulses or warning	channel 1 (M) warning (W) channel 2 (K)
Size	$\frac{1}{4}$ -unit 19"-system
Height	178 mm
Weight	about 2500 grammes

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Anticoincidence-coincidence-amplifier unit, type AKV 2

Description

The logical circuit type AKV 2 (DP) is used for the measurement of anticoincident and coincident pulses, generated in our tritium detectors of the series DDH; with it may be carried out, beside of the specific measurement of tritium, also the specific measurement of higher-energy isotopes in air. The circuitry may be modified according as the problem, whereby the designation alters. If merely the tritium content of air is to be measured or only the need for a fair anticoincidence measurement is present, as with our large area counters of the series DDA, so the unit bears the designation AV 2. In the same unit also one-channel amplifiers, type V 12, two-channel amplifiers, type V 22, and an anticoincidence circuitry, type GV 2, effecting the two channels, are built in. In any case the panel shows the block diagram of the built-in circuitry. The respective detectors or their systems are

connected to the high voltage output on the panel. The high voltage is fed in a high voltage input on the rear from the unit, type HSP 2; the power supply follows from the unit, type STW 2. By request the power supply can be built in the unit for special cases.

The amplifiers of this unit guarantee due to inverse feedback high constancy and steadiness of the overriding. The latter is important, since very high amplitudes get to the amplifier input, when both alpha and beta radiations are measured.

Each circuitry is built in a $\frac{1}{4}$ -unit in 19"-system. The unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. All parts are easily accessible and interchangeable.

Technical data

(subject to change without notice)

Input sensitivity of a channel	0,001 V
Input resistance	20 pF/500 000 ohms
Output voltage of a channel or a logical circuitry	5 V pp, negativ
Standard output pulses, duration	2 usec.
Output resistance	50 ohms
Integrated amplifier with inverse feedback	
Fail-safe for pulses to	5 V negative
Power supply from	STW 2
Size	$\frac{1}{4}$ -unit, 19"-system
Height	178 mm
Weight (AKV 2)	about 1500 grammes

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
Ratemeter unit, type ERM 2

Application

The ratemeter unit is intended for an accurate analogous indication of pulse rates, that are fed in it through a preamplifying circuit (e. g. type AV 2). Since it may be operated, beside of operating by hand, in addition by full-automatic operating or remote control, it also is suited for monitoring and controls free from maintenance.

Construction

The device is built in a $\frac{1}{2}$ -unit in 19"-system. The unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. All parts are easily accessible and interchangeable. The controls on the panel are: electric zero shift, switch for connection to automatic range switching, 8 push buttons for measuring ranges, 6 push buttons for time constants, two adjustable warning thresholds on the measuring contactor

Performance

(see block diagramm)

The pulse rates fed in the ratemeter give a DC voltage on the capacitor C, that determines the time constants. This DC voltage is amplified by an indication-record amplifier and passed on to the record output on the one hand and to the measuring contactor on the other hand. The latter indicates the measuring value and has minimum and maximum contacts for the adjustment of warning thresholds. An electric zero shift permits to suppress possible deflections of the indicating measuring contactor, e. g. originating from disturbing radiations or background. Nevertheless the accuracy of indication and the sensitivity of the selected measuring range are kept. Thus a compensation is possible at constant level of disturbing radiations.

Technical data

(subject to change without notice)

Input pulses from	$\frac{1}{4}$ -unit V 12, AV 2, AKV 2
8 linear measuring ranges	100...300, 1000...3000, 10000...30000, 100000...300000 pulses per minute
Accuracy of indication	2 %
Deviation from linearity within a measuring range	1 %
6 time constants, adjustable by push buttons	2, 4, 16, 32, 64, 128 sec.
Indicator with 2 independent contacts for adjustment of a maximum and minimum warning threshold	measuring contactor
Modes of operation: Measuring ranges to select by hand, remote control, and automatic range switching	
Electric zero shift, disconnectable	— 100 %
Possibility of linear overriding within a range at zero suppression	100 %
Record output	max. 10 V, max. 0,01 A against ground
Power supply from	$\frac{1}{4}$ -unit STW 2
Size	$\frac{1}{2}$ -unit, 19"-system
Height	178 mm
Weight	about 4300 grammes


Block diagramm ERM 2

(from top left to bottom right)

ratemeter measuring channel — time constant C — indication-record amplifier — range switching — measuring contactor — minimum-maximum threshold

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
Automatic range switch, type ABU 2

Application

The automatic range switch (DPA.) may be used in connection with the difference ratemeter, type DRM 2, and the ratemeter, type ERM 2. The automatic control permits an operation-free switching and remote control of the ratemeters. Hence the automatic range switch is brought to bear in continuous measurements of air and water, in measurement devices, separated relating to space, and control rooms. Also it permits the release of control processes.

Construction

The device is built in a $\frac{1}{4}$ -unit in 19"-system. The unit is enclosed and all voltages are connected through rack and panel connectors on the rear of the unit. The circuitry is largely carried out with plug boards, printed-circuit boards, and IC's. All parts are easily accessible and interchangeable.

Performance

(see block diagramm)

The indication-recorder-output voltage, generated in the ratemeter, is fed in a min/max-switch. This switch actuates, according as the account of the fed voltage, through a logical circuitry a forward-backward counter, that counts forwards or backwards according as the position of the min/max-switch. The switching speed is determined by a time pulse generator. A decoder, attached to the forward-backward counter, operates the range relays through current amplifiers. The range relays have outputs for the measuring ranges of the ratemeter, warning, and measuring range remote indication. For example, if the output voltage of the ratemeter is too high, the ratemeter is switched in the next after range.

The time pulse generator provides for switching as long till the output voltage is again between maximum and minimum of the values positioned on the switch.

Technical data

(subject to change without notice)

Automatic switching range for 8 measuring ranges

Max/min-switch, continuously variable for voltages 0...10 V

Minimum-maximum contacts adjustable for any range

Range relay comb relay

Built-in voltage step function for range remote indication

Power supply from STW 2
 ± 15 V and ± 5 V against ground

Size $\frac{1}{4}$ -unit, 19"-system

Height 178 mm

Weight about 2000 grammes

Block diagramm ABU 2

(from top left to bottom right)

minimum switch — logical circuitry — forward-backward counter — decoder — ratemeter output voltage — maximum switch — time pulse generator — range relays — outputs for: measuring ranges, warning, range remote indication